



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Tube,sterile,screw cap,5.0 ml,BOX/500

**Product Code :**

LBNY-0007-22600061



**Description :**

Applicable for molecular/PCR & sterile microbiological works. Various clinical laboratory usage conditions including; Sample cold storage, sample preparation/extractions with water bath, incubator and centrifuge.

**Technical Specification :**

Polypropylene material  
5.0ml capacity  
Sterile  
Free of DNA, RNase, DNase & PCR inhibitors  
Conical base  
Screw cap  
Graduated  
Withstands centrifugation force up to 25,000 x g.

**Weight & volume:**

Estimated volume: 42.00dm<sup>3</sup>

Estimated weight: 4.49kg.

---

# Laboratory Instrument India

