



Email : sales@laboratoryinstrumentindia.com

Product Name :

Tube,sterile,screw cap,5.0 ml,BOX/500

Product Code :

LBNY-0007-22600061



Description :

Applicable for molecular/PCR & sterile microbiological works. Various clinical laboratory usage conditions including; Sample cold storage, sample preparation/extractions with water bath, incubator and centrifuge.

Technical Specification :

Polypropylene material
5.0ml capacity
Sterile
Free of DNA, RNase, DNase & PCR inhibitors
Conical base
Screw cap
Graduated
Withstands centrifugation force up to 25,000 x g.

Weight & volume:

Estimated volume: 42.00dm³

Estimated weight: 4.49kg.

Laboratory Instrument India

