



Email : sales@laboratoryinstrumentindia.com

Product Name :
Tube, snapcap, conic, 1.5ml, BOX-1000

Product Code :
LBNY-0007-22600062



Description :

Various clinical laboratory uses including Centrifugation and short-term sample storage.

Technical Specification :

Polypropylene material
1.5ml capacity
Secure snap cap attached to the tube
Free of DNA, RNase, DNase & PCR inhibitors
Conical base
Non-sterile
Graduated
Compatible with the standard microcentrifuges
Withstands centrifugation force up to 25,000 x g.

Weight & volume:
Estimated volume: 11.57dm³
Estimated weight: 1.35kg.

Laboratory Instrument India

