



Email : sales@laboratoryinstrumentindia.com

Product Name :

Tube, snapcap, sterile, conic, 1.5ml, BOX-200

Product Code :

LBNY-0007-22600063

**Description :**

Applicable for molecular/PCR & sterile microbiological works. Various clinical laboratory uses including Centrifugation and short-term sample storage.

Technical Specification :

Polypropylene material
1.5ml capacity
Secure snap cap attached to the tube
Withstands centrifugation force up to 25,000 x g
Sterile
Free of DNA, RNase, DNase & PCR inhibitors
Conical base
Graduated
Compatible with the standard microcentrifuges.

Weight & volume:

Estimated volume: 2.85dm³

Estimated weight: 0.3kg.

Laboratory Instrument India

