



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Tube, snapcap, sterile, conic, 1.5ml, BOX-200

**Product Code :**

LBNY-0007-22600063



**Description :**

Applicable for molecular/PCR & sterile microbiological works. Various clinical laboratory uses including Centrifugation and short-term sample storage.

**Technical Specification :**

Polypropylene material  
1.5ml capacity  
Secure snap cap attached to the tube  
Withstands centrifugation force up to 25,000 x g  
Sterile  
Free of DNA, RNase, DNase & PCR inhibitors  
Conical base  
Graduated  
Compatible with the standard microcentrifuges.

**Weight & volume:**

Estimated volume: 2.85dm<sup>3</sup>

Estimated weight: 0.3kg.

---

# Laboratory Instrument India

