



Email : sales@laboratoryinstrumentindia.com

Product Name :
Centrifuge, haematocrit

Product Code :
LBNY-0001-2160001



Description :

Intended for clinical laboratory use; Haematocrit/PCV, Blood separation etc.

Technical Specification :

Automatic brake and lid interlock
Automated speed approx. up to 15000 rpm
Timer 1-99 minutes
Backlit LCD display, digital operable by touch /press allowing responsive operation – even when wearing gloves
Capable of handling capillary tubes.
Power requirements 110-220V/50-60Hz, with voltage surge protection
High sample throughput, multi-safety control, personal program function
Plastic housing, metal lid, deep metal chamber which securely holds the rotor
Rotor holds 24 capillary tubes of 75mm.

Weight & volume:
Estimated volume:0.023m3
Estimated weight: 10kg.

Laboratory Instrument India

