



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Centrifuge, haematocrit

**Product Code :**  
LBNY-0001-2160001



**Description :**

Intended for clinical laboratory use; Haematocrit/PCV, Blood separation etc.

**Technical Specification :**

Automatic brake and lid interlock  
Automated speed approx. up to 15000 rpm  
Timer 1-99 minutes  
Backlit LCD display, digital operable by touch /press allowing responsive operation – even when wearing gloves  
Capable of handling capillary tubes.  
Power requirements 110-220V/50-60Hz, with voltage surge protection  
High sample throughput, multi-safety control, personal program function  
Plastic housing, metal lid, deep metal chamber which securely holds the rotor  
Rotor holds 24 capillary tubes of 75mm.

Weight & volume:  
Estimated volume:0.023m3  
Estimated weight: 10kg.

---

# Laboratory Instrument India

