



Email : sales@laboratoryinstrumentindia.com

Product Name :
Centrifuge, basic serology, 24x5ml

Product Code :
LBNY-0001-2160002



Description :

Blood sample separation for clinical laboratory tests (serum, plasma, etc).

Technical Specification :

Bench top type
Metal housing, aerosol-proof-lid, stainless steel chamber.
Automatic imbalance recognition.
Large LCD screen with clear indication of running condition including time and speed
Lid locking and holding, emergency lid lock release
Power requirements 220V / 50Hz
Rotor speed up to 6000rpm
Timer 1-99 minutes.

Weight & volume:
Estimated volume: 0.042m³
Estimated weight 23.00kg.



Laboratory Instrument India

