



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Centrifuge, basic serology, 24x5ml

**Product Code :**  
LBNY-0001-2160002



**Description :**

Blood sample separation for clinical laboratory tests (serum, plasma, etc).

**Technical Specification :**

Bench top type  
Metal housing, aerosol-proof-lid, stainless steel chamber.  
Automatic imbalance recognition.  
Large LCD screen with clear indication of running condition including time and speed  
Lid locking and holding, emergency lid lock release  
Power requirements 220V / 50Hz  
Rotor speed up to 6000rpm  
Timer 1-99 minutes.

Weight & volume:  
Estimated volume: 0.042m<sup>3</sup>  
Estimated weight 23.00kg.



**Laboratory Instrument India**

