



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Centrifuge, complete w/accessories

**Product Code :**  
LBNY-0001-2160003



**Description :**

Suitable for use in resource-limited settings. For the separation of heterogeneous mixtures or samples in clinical laboratory applications.

**Technical Specification :**

Bench top type  
Large LCD screen with clear indication of running condition including time and speed  
Metal housing, aerosol-proof-lid, stainless steel chamber.  
Automatic imbalance recognition.  
Rotor speed up to 6000rpm  
Timer 1-99 minutes  
Power requirements 220V / 50Hz  
Lid locking and holding, emergency lid lock release

Weight & volume:  
Estimated volume: 0.042m<sup>3</sup>  
Estimated weight 23.00kg.

---

# Laboratory Instrument India

