



Email : sales@laboratoryinstrumentindia.com

Product Name :
Centrifuge, complete w/accessories

Product Code :
LBNY-0001-2160003



Description :

Suitable for use in resource-limited settings. For the separation of heterogeneous mixtures or samples in clinical laboratory applications.

Technical Specification :

Bench top type
Large LCD screen with clear indication of running condition including time and speed
Metal housing, aerosol-proof-lid, stainless steel chamber.
Automatic imbalance recognition.
Rotor speed up to 6000rpm
Timer 1-99 minutes
Power requirements 220V / 50Hz
Lid locking and holding, emergency lid lock release

Weight & volume:
Estimated volume: 0.042m³
Estimated weight 23.00kg.

Laboratory Instrument India

