



Email : sales@laboratoryinstrumentindia.com

Product Name :
Centrifuge, manual, 4 positions

Product Code :
LBNY-0001-2160004



Description :

Used to centrifuge clinical samples at field facilities where there is no electrical power with a clamping device for simple table mounting.

Technical Specification :

Bench clamp with adaptable screw
Four plastic buckets for test tubes
Manual, hand operated
Metal base
Self-lubrication system
Shaft speed up to 1298RCF or 3000 rpm
Metal dismountable gear support structure
Swing-out rotator, 4-places at 90 degrees.

Weight & volume:
Estimated volume: 0.008m³
Estimated weight: 1.00kg.

Laboratory Instrument India

