



Email : sales@laboratoryinstrumentindia.com

Product Name :
Shaker, orbital

Product Code :
LBNY-0001-2160005



Description :

Mixing of reagents in bottles, flasks and beakers in horizontal orbital motion.

Technical Specification :

Microprocessor controlled shaking platform providing smooth and quiet horizontal orbital motion for mixing in bottles, flasks, and beakers

Continuous or timed operation with automatic switch-off

Power requirements 110-220V/50Hz-60Hz, with voltage surge protection.

Automatic load balancing system

Audible signal accompanies automatic switch-off to indicate the completion of a set time period

Adjustable shaking speed 50 to 450 rpm

Amplitude 3mm

Orbit up to 10mm

Platform size: 10x15cm, accommodates flasks, beakers, dishes and Petri-dishes

Load up to 4.5kg.

Weight & volume:

Estimated volume: 0.005m³

Estimated weight: 3.00kg.

Laboratory Instrument India

