



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Rotator, blood specimen

**Product Code :**  
LBNY-0001-2160006



**Description :**

Rotator mixer for blood samples to ensure a homogenous mixture for clinical laboratory applications.

**Technical Specification :**

Rocking motion and rotation around a horizontal axis  
Designed for continuous or timed quiet operation  
Adjustable speed from 5 to 60 rpm  
Capacity 6 rollers with 22 blood tubes of 10 -16mm diameter (1.5 ml-15ml tubes)  
Built-in digital timer mode.

Weight & volume:  
Estimated volume:0.026m<sup>3</sup>  
Estimated weight: 7.00kg.



**Laboratory Instrument India**