



Email : sales@laboratoryinstrumentindia.com

Product Name :
Vortex, test tube

Product Code :
LBNY-0001-2160007



Description :

For clinical laboratory applications. Used to ensure thorough and vigorous mixing of samples and reagents.

Technical Specification :

Compact design mixer/vibrator of a single tube
Continuous adjustable speed up to 3300 rpm.
Circular horizontal orbit up to 5mm.
Operates continuously or starts mixing when rubber top is pressed down with tube
Power requirements 110-220V / 50Hz-60Hz, with voltage surge protection.
Rubber top cavity for tubes.
Non-skid feet prevent the vortex mixer from sliding.

Weight & volume:
Estimated volume: 0.003m³
Estimated weight: 3.5kg.



Laboratory Instrument India

