



Email : sales@laboratoryinstrumentindia.com

Product Name :
Water bath, 7 L

Product Code :
LBNY-0001-2160009



Description :

For clinical laboratory applications. Equipment, bench-top, filled with water, and used as a controlled temperature environment for incubating samples and/or tests.

Technical Specification :

Minimum water carrying capacity 7 liters
Inbuilt overheating protection system
Power requirements 230V with voltage surge protection.
Temperature range from 25C to 100C
Variations within the bath approx. 0.1C
Stainless steel bath interior, cover, sieve plate, and heating element.
Visual and acoustic alarm for low-level and over/under temperature.

Weight & volume:
Estimated volume: 0133m3
Estimated weight: 11.00kg.



Laboratory Instrument India

