



Email : sales@laboratoryinstrumentindia.com

Product Name :
Distillation unit, 2 L/h, with tank

Product Code :
LBNY-0001-21600010



Description :

Double-walled housing. Heating element, evaporator, storage tank, and condenser made of stainless steel
Distilled water to comply with the International Pharmacopoeia, free from pyrogen and have a conductivity of 2 μ S/cm at 25oC
Automatic system with minimum production 2 L/h
Suitable for both bench and wall mounting
Operates from raw potable feed water, 3 - 7 bar, 4oC to 35oC
Built-in storage tank capacity approx. 4L minimum
Water saving through automatic water cut-off when storage tank is full
Thermostatic overheat protection
Mechanism to automatically restart the distiller when distilled water is withdrawn
Electronic contamination monitoring of the still with distillation stop and visual cleaning instruction
Thermostatic low water cut off protects the heating element
Solenoid valve to control the water supply
Power requirements 230V / 50-60Hz.

Technical Specification :

Weight & volume:
Estimated volume: 0.15m³
Estimated weight: 24.00kg.



Laboratory Instrument India