



Email : sales@laboratoryinstrumentindia.com

Product Name :
PCR Plates Centrifuge

Product Code :
LBNY-0001-21600016



Description :

For quick and easy centrifugation of samples in PCR plates, strips, tubes or microplates.

Technical Specification :

Capacity: Micro-plates 384-PCR/96-PCR/Elisa plates x 8 pcs; Deep well plates x 4 pcs
Timer can be set up to 99 min. for short run/ continuous run
Pre-programmed and user-defined programmes
LCD display
Force: 2200 x G
Speed 300-6000rpm. Increment:10rpm
Power requirements 220V / 50Hz, with voltage surge protection.



Laboratory Instrument India