



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
PCR Plates Centrifuge

**Product Code :**  
LBNY-0001-21600016



**Description :**

For quick and easy centrifugation of samples in PCR plates, strips, tubes or microplates.

**Technical Specification :**

Capacity: Micro-plates 384-PCR/96-PCR/Elisa plates x 8 pcs; Deep well plates x 4 pcs  
Timer can be set up to 99 min. for short run/ continuous run  
Pre-programmed and user-defined programmes  
LCD display  
Force: 2200 x G  
Speed 300-6000rpm. Increment:10rpm  
Power requirements 220V / 50Hz, with voltage surge protection.



**Laboratory Instrument India**