



Email : sales@laboratoryinstrumentindia.com

Product Name :
ICE Maker

Product Code :
LBNY-0001-21600017



Description :

Daily Ice Making Capacity: 20kg/24h
Water Tank Capacity: 1L
Automatic Overflow Protection
Ice Shape: cube of approx. 13gm each
Installation: Freestanding Air Cooled
Minimum Cycle Time for First Ice Made: 30 minutes
Nonremovable ice bin
Voltage: 220V-230V/50Hz
Ice Bucket Capacity: 4kg (approx. 310 ice cubes).

Technical Specification :

Weight & volume:
Estimated volume: 0.24m³
Estimated weight: 35.00kg.



Laboratory Instrument India

