



Email : sales@laboratoryinstrumentindia.com

Product Name :
Micro Centrifuge, 24x2ml

Product Code :
LBNY-0001-21600018



Description :

Intended to provide centrifugal separation in 2.0ml microtubes in the different fields of medical laboratories.

Technical Specification :

Fixed angle rotor with capacity for 24 standard 2.0ml microcentrifuge tubes/cryovials
Backlit LCD display, digitally operable by touch /press allowing responsive operation – even when wearing gloves
Automated speed Max RCF up to 25000 x g
Stainless steel chamber
Power requirements 200-240V/50Hz
Timer 1 - 99 minutes
Automatic brake and lid interlock
Inbuilt temperature control. Adjustable refrigerated centrifugation setting from -20oC to 40oC.



Laboratory Instrument India