



Email : sales@laboratoryinstrumentindia.com

Product Name :
Pipette, digital, 2-20 ul

Product Code :
LBNY-0002-2160001



Description :

Used for multiple techniques to accurately dispense required volumes of liquid.
Single channel, adjustable pipetting volumes 2-20ul
Autoclavable, and Resistant to UV sterilization
Sturdy and chemical resistant parts.
A handle with digital setting and read-out of delivered volumes
Built-in tips ejector
Manual air displacement micropipette
Compatible with universal fit pipette tips ranging 2-20ul.

Technical Specification :

Weight & volume:
Estimated Volume: 2.20dm³
Estimated Weight: 0.105kg.



Laboratory Instrument India