



Email : sales@laboratoryinstrumentindia.com

Product Name :
Pipette, digital, 10-100 ul

Product Code :
LBNY-0002-2160002



Description :

Used for multiple techniques to accurately dispense required volumes of liquid
Single channel, adjustable pipetting volumes 10-100ul
Built-in tips ejector
Compatible with universal fit pipette tips ranging 10-100ul
Autoclavable, and Resistant to UV sterilization
Manual air displacement micropipette
A handle with digital setting and read-out of delivered volumes
Sturdy and chemical resistant parts.

Technical Specification :

Weight & volume:
Estimated Volume: 2.20dm³
Estimated Weight: 0.106kg.



Laboratory Instrument India