



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Pipette, digital, 100-1000 ul

**Product Code :**  
LBNY-0002-2160004



**Description :**

Used for multiple techniques to accurately dispense required volumes of liquid.  
A handle with digital setting and read-out of delivered volumes  
Built-in tips ejector  
Manual air displacement micropipette  
Compatible with universal fit pipette tips ranging 100-1000ul  
Autoclavable, and Resistant to UV sterilization  
Single channel, adjustable pipetting volumes 100-1000ul  
Sturdy and chemical resistant parts.

**Technical Specification :**

Weight & volume:  
Estimated Volume: 2.20dm<sup>3</sup>  
Estimated Weight: 0.123kg.



**Laboratory Instrument India**