



Email : sales@laboratoryinstrumentindia.com

Product Name :
Pipette, digital, 8 channel, 20-200ul

Product Code :
LBNY-0002-2160006



Description :

Used for multiple techniques to accurately dispense required volumes of liquid.
8-channel, adjustable pipetting volumes 20-200ul
Compatible with universal fit pipette tips ranging 20-200ul
Autoclavable, and Resistant to UV sterilization
Manual air displacement micropipette
A handle with digital setting and read-out of delivered volumes
Built-in tips ejector
Sturdy and chemical resistant parts.

Technical Specification :

Weight & volume:
Estimated Volume: 2.20dm³
Estimated Weight: 0.209kg.



Laboratory Instrument India