



Email : sales@laboratoryinstrumentindia.com

Product Name :

Pipette, digital, 8 channel, 100-1200ul

Product Code :

LBNY-0002-2160007



Description :

Used for multiple techniques to accurately dispense required volumes of liquid
8-channel, adjustable pipetting volumes 100-1200ul
Compatible with universal fit pipette tips ranging 100-1000ul
Autoclavable, and Resistant to UV sterilization
Manual air displacement micropipette
A handle with digital setting and read-out of delivered volumes
Built-in tips ejector
Sturdy and chemical resistant parts.

Technical Specification :

Weight & volume:
Estimated Volume: 0.003m³
Estimated Weight: 0.258kg.



Laboratory Instrument India

