



Email : sales@laboratoryinstrumentindia.com

Product Name :

Pipette, digital, 12 channels, 5-50 ul

Product Code :

LBNY-0002-2160008



Description :

Used for multiple techniques to accurately dispense required volumes of liquid
12-channel, adjustable pipetting volumes 5-50ul
Compatible with universal fit pipette tips ranging 5-50ul
Autoclavable, and Resistant to UV sterilization
A handle with digital setting and read-out of delivered volumes
Manual air displacement micropipette
Built-in tips ejector
Sturdy and chemical resistant parts.

Technical Specification :

Weight & volume:
Estimated Volume: 2.20dm³
Estimated Weight: 0.243kg.



Laboratory Instrument India

