



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Pipette, digital, 12 channels, 5-50 ul

**Product Code :**  
LBNY-0002-2160008



**Description :**

Used for multiple techniques to accurately dispense required volumes of liquid  
12-channel, adjustable pipetting volumes 5-50ul  
Compatible with universal fit pipette tips ranging 5-50ul  
Autoclavable, and Resistant to UV sterilization  
A handle with digital setting and read-out of delivered volumes  
Manual air displacement micropipette  
Built-in tips ejector  
Sturdy and chemical resistant parts.

**Technical Specification :**

Weight & volume:  
Estimated Volume: 2.20dm<sup>3</sup>  
Estimated Weight: 0.243kg.



**Laboratory Instrument India**

