



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Pipette, digital, 12 channels, 20-200 ul

**Product Code :**

LBNY-0002-2160009



**Description :**

Used for multiple techniques to accurately dispense required volumes of liquid  
12-channel, adjustable pipetting volumes 20-200ul  
Compatible with universal fit pipette tips ranging 20-200ul  
Autoclavable, and Resistant to UV sterilization  
Manual air displacement micropipette  
A handle with digital setting and read-out of delivered volumes  
Sturdy and chemical resistant parts  
Built-in tips ejector.

**Technical Specification :**

Weight & volume:  
Estimated Volume: 2.20dm<sup>3</sup>  
Estimated Weight: 0.248kg.



**Laboratory Instrument India**

