



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Pipette, digital, 12 channel, 50-1200ul

**Product Code :**

LBNY-0002-21600010



**Description :**

Used for multiple techniques to accurately dispense required volumes of liquid  
Compatible with universal fit pipette tips ranging 100-1000µl  
Autoclavable, and Resistant to UV sterilization  
12-channel, adjustable pipetting volumes 50-1200µl  
A handle with digital setting and read-out of delivered volumes  
Manual air displacement micropipette  
Built-in tips ejector  
Sturdy and chemical resistant parts.

**Technical Specification :**

Weight & volume:  
Estimated Volume: 0.003m<sup>3</sup>  
Estimated Weight: 0.365kg.



**Laboratory Instrument India**

