



Email : sales@laboratoryinstrumentindia.com

Product Name :

Pipette, digital, 12 channel, 50-1200ul

Product Code :

LBNY-0002-21600010



Description :

Used for multiple techniques to accurately dispense required volumes of liquid
Compatible with universal fit pipette tips ranging 100-1000 μ l
Autoclavable, and Resistant to UV sterilization
12-channel, adjustable pipetting volumes 50-1200 μ l
A handle with digital setting and read-out of delivered volumes
Manual air displacement micropipette
Built-in tips ejector
Sturdy and chemical resistant parts.

Technical Specification :

Weight & volume:
Estimated Volume: 0.003m³
Estimated Weight: 0.365kg.



Laboratory Instrument India

