



Email : sales@laboratoryinstrumentindia.com

Product Name :

Box, specimen transp., 2L and 4L, set/2

Product Code :

LBNY-0001-2170001



Description :

For transporting clinical samples and products that you want to keep clean and dry under testing environmental conditions.

Technical Specification :

Can hold samples and tubes, for transportation to the clinical laboratory
2 L and 4 L.

Snap lock securely clasps the lid.

Sturdy carry handle integrated with lid.

Chemical resistant and autoclavable plastic.

Set of 1 box supplied with tube rack sizes are specified below:

1 x rack for 40 standard test tubes or 10 specimen containers.

1 x rack for 20 standard test tubes.

Weight – Volume:

Estimated volume: 0.04m³.

Estimated weight: 2 kg.

Laboratory Instrument India

