



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Tourniquet, with buckle

**Product Code :**  
LBNY-0001-2170002



**Description :**

To help locate and define peripheral veins to achieve successful and safe venepuncture.  
Self-applicable with a quick-release system.  
Reusable, washable.  
With easy snap-lock/press-release buckle, made of ABS plastic or equivalent durable material.  
Made of durable elastic belt material.  
Color: Green.  
Length: 40 cm.  
Width: 3 cm.

**Technical Specification :**

Weight & volume:  
Estimated Volume: 0.0005m<sup>3</sup>.  
Estimated Weight: 0.15 Kg.



**Laboratory Instrument India**

