



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Paper, dry blood spot, box/100

**Product Code :**  
LBNY-0001-2170003



**Description :**

Designed for blood sample collection, preserves nucleic acid in the dried blood spots for subsequent laboratory tests.

To collect and transport dry blood specimens for early infant diagnosis of HIV

Dry blood spot (DBS) absorbent filter paper.

Stability: Minimum 24 months for DNA.

DNA, DNase, RNase, & PCR inhibitor.

With 5 distinct pre-marked spots.

Spots are pre-perforated / pre-punched, facilitating processing in PCR applications.

**Technical Specification :**

Weight & volume:

Estimated Volume: 0.0005m<sup>3</sup>.

Estimated Weight: 0.15 Kg.



**Laboratory Instrument India**

