




Email : sales@laboratoryinstrumentindia.com

Product Name :
Rack, drying DBS cards, 10 positions

Product Code :
LBNY-0001-2170004



Description :

To hold dry blood spot cards during air drying for early infant diagnosis of HIV.
Provides horizontal position for 10 cards.
Sturdy base provides optimal stability.
Self-supporting vertical rack for drying blood spot filter cards.
Material: Cardboard.

Technical Specification :

Weight & volume:
Estimated Volume: 0.0006m3.
Estimated Weight: 0.04 kg.



Laboratory Instrument India