



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Tube, capillary, heparin, box/1000

**Product Code :**  
LBNY-0001-2170005



**Description :**

Used for multiple tests.  
Heparin coated internally  
Capillary glass tube  
Length approx. 75 mm  
Diameter approx. 1.2 mm  
Seal-packing of smaller units (such as 100 tubes).

**Technical Specification :**

Weight & volume:  
Estimated Volume: 0.009m3.  
Estimated Weight: 1 Kg.



**Laboratory Instrument India**