



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Tube, capillary, EDTA, box/1000

**Product Code :**  
LBNY-0001-2170006



**Description :**

Used for multiple tests.  
EDTA coated internally.  
Capillary glass tube.  
Length approx. 75 mm.  
Diameter approx. 1.2 mm.

**Technical Specification :**

Shelf life:  
48 months.

Weight & volume:  
Estimated Volume: 0.009m3.  
Estimated Weight: 1 Kg.



**Laboratory Instrument India**

