



Email : sales@laboratoryinstrumentindia.com

Product Name :

Tube, vacuum, EDTA, 2 ml, box/100

Product Code :

LBNY-0001-2170007



Description :

Blood collection for clinical laboratory tests.
Blood collection vacuum tube.
Material: clear non-fragile PET plastic and leak-proof rubber cap.
Plastic shielded rubber cap to protect personnel from direct contact with blood.
Capped, with needle pierceable vacuum seal.
Additive EDTA anticoagulant
Standard color-coded cap.
Applicable for automated analysers.
Sterile.
Blood draw volume 2.0ml.
Pen/Marker writable label space.

Technical Specification :

Shelf life:
16 months.

Weight & volume:
Estimated Volume: 0.002m3.
Estimated Weight: 0.4 Kg.



Laboratory Instrument India