



Email : sales@laboratoryinstrumentindia.com

Product Name :

Tube, vacuum, EDTA, 4 ml, box/100

Product Code :

LBNY-0001-2170008



Description :

Blood collection for clinical laboratory tests.
Capped, with a needle-pierceable vacuum seal.
Blood collection vacuum tube.
Material: clear nonfragile PET plastic and leak-proof rubber cap.
Plastic shielded rubber cap to protect personnel from direct contact with blood.
Applicable for automated analyzers.
Pen/Marker writable label space.
Blood draw volume 4.0ml.
Additive EDTA anticoagulant.
Standard color-coded cap.

Technical Specification :

Shelf life:
16 Months.

Weight and Volume:
Estimated volume: 0.0034m3.
Estimated weight: 0.7 Kg.

Laboratory Instrument India

