



Email : sales@laboratoryinstrumentindia.com

Product Name :

Tube, vacuum, EDTA, 6 ml, box/100

Product Code :

LBNY-0001-2170009



Description :

Blood collection vacuum tube.

Material: clear nonfragile PET plastic and leak-proof rubber cap.

Plastic shielded rubber cap to protect personnel from direct contact with blood.

Sterile.

Blood draw volume 6.0ml.

Capped, with needle-pierceable vacuum seal cap.

Applicable for automated analyzers.

Pen/Marker writable label space.

Additive EDTA anticoagulant.

Standard color-coded cap.

Technical Specification :

Shelf life:

18 months.

Weight and Volume:

Estimated Volume: 0.00495m³.

Estimated Weight: 0.825 kg.

Laboratory Instrument India

