



Email : sales@laboratoryinstrumentindia.com

Product Name :

Tube, vacuum, serum gel, 4 ml, box/100

Product Code :

LBNY-0001-21700010



www.laboratoryinstrumentindia.com

Description :

Serum collection for clinical laboratory tests.

Blood draw volume 3.5ml.

Material: clear nonfragile PET plastic and leak-proof rubber cap.

Plastic shielded rubber cap to protect personnel from direct contact with blood.

Blood collection vacuum tube.

Capped, with needle-pierceable vacuum seal cap.

Additive; Serum separation gel with clot activator.

Pen/Marker writable label space.

Standard color-coded cap.

Technical Specification :

Shelf life:

15 months.

Weight and Volume:

Estimated volume: 0.0034m³.

Estimated weight: 0.77kg.



Laboratory Instrument India