




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Tube, vacuum, serum gel, 6 ml, box/100

**Product Code :**

LBNY-0001-21700011



**Description :**

Blood collection vacuum tube.  
Blood draw volume 6.0ml.  
Material: clear nonfragile PET plastic and leak-proof rubber cap.  
Plastic shielded rubber cap to protect personnel from direct contact with blood.  
Capped, with needle-pierceable vacuum seal cap.  
Standard color-coded cap.  
Pen/Marker writable label space.  
Additive; Serum separation gel with clot activator.

**Technical Specification :**

Shelf life:  
16 months.

Weight and volume:  
Estimated volume:0.004954 m3.  
Estimated Weight: 0.9kg.

---

# Laboratory Instrument India

