




Email : sales@laboratoryinstrumentindia.com

Product Name :

Tube, vacuum, serum gel, 6 ml, box/100

Product Code :

LBNY-0001-21700011



Description :

Blood collection vacuum tube.
Blood draw volume 6.0ml.
Material: clear nonfragile PET plastic and leak-proof rubber cap.
Plastic shielded rubber cap to protect personnel from direct contact with blood.
Capped, with needle-pierceable vacuum seal cap.
Standard color-coded cap.
Pen/Marker writable label space.
Additive; Serum separation gel with clot activator.

Technical Specification :

Shelf life:
16 months.

Weight and volume:
Estimated volume:0.004954 m3.
Estimated Weight: 0.9kg.

Laboratory Instrument India

