



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Tube, vacuum, plain/dry, 4 ml, box/100

**Product Code :**

LBNY-0001-21700012



**Description :**

Serum collection for clinical laboratory tests

Blood collection vacuum tube

Blood draw volume 4.0ml

Material: clear nonfragile PET plastic and leak-proof rubber cap

Plastic shielded rubber cap to protect personnel from direct contact with blood

Capped, with needle pierceable vacuum seal cap

Sterile

Plain/dry, No additive

Pen/Marker writable label space.

**Technical Specification :**

Shelf life:

16 months.

Weight – Volume:

Estimated volume: 0.00344 m3.

Estimated weight: 0.9 kg.

---

# Laboratory Instrument India

