



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Tube, vacuum, plain/dry, 6 ml, box/100

**Product Code :**

LBNY-0001-21700013



**Description :**

Serum collection for clinical laboratory tests.

Blood collection vacuum tube.

Blood draw volume 6.0ml.

Material: clear nonfragile PET plastic and leak-proof rubber cap.

Plastic shielded rubber cap to protect personnel from direct contact with blood.

Capped, with needle-pierceable vacuum seal cap.

Plain/dry, No additive.

Pen/Marker writable label space.

Standard color-coded cap.

**Technical Specification :**

Shelf life:

18 months.

Weight and Volume:

Estimated volume: 0.00495 m3.

Estimated weight: 0.825 kg.

---

# Laboratory Instrument India

