



Email : sales@laboratoryinstrumentindia.com

Product Name :
Needle, vacuum tube, 20 G, box/100

Product Code :
LBNY-0001-21700014



Description :

Blood collection needle for use with vacuum tubes. For use in clinical phlebotomy.
Blood collection needle.
Sterile, individually packed in a twist-apart plastic container.
Material; Medical grade stainless steel.
Length approx. 40 mm
Diameter 20 G (0.9 mm)
Applicable for multi-sample draw.
Needle fits (screw in and out) needle holder.
Pointed at both ends, with rubber sleeve on the posterior end needle.

Technical Specification :

Shelf life:
61 months.

Weight and Volume:
Estimated volume: 0.0011m3.
Estimated weight: 0.26 kg.



Laboratory Instrument India