



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Needle, vacuum tube, 20 G, box/100

**Product Code :**

LBNY-0001-21700014



**Description :**

Blood collection needle for use with vacuum tubes. For use in clinical phlebotomy.

Blood collection needle.

Sterile, individually packed in a twist-apart plastic container.

Material; Medical grade stainless steel.

Length approx. 40 mm

Diameter 20 G (0.9 mm)

Applicable for multi-sample draw.

Needle fits (screw in and out) needle holder.

Pointed at both ends, with rubber sleeve on the posterior end needle.

**Technical Specification :**

Shelf life:

61 months.

Weight and Volume:

Estimated volume: 0.0011m<sup>3</sup>.

Estimated weight: 0.26 kg.



**Laboratory Instrument India**