



Email : sales@laboratoryinstrumentindia.com

Product Name :

Needle, vacuum tube, 22 G, box/100

Product Code :

LBNY-0001-21700015



Description :

Applicable for multi-sample draw
Needle fits (screw in and out) needle holder
Sterile, individually packed in a twist-apart plastic container
Blood collection needle for use with vacuum tubes. For use in clinical phlebotomy.
Material; Medical grade stainless steel
Blood collection needle
Pointed at both ends, with rubber sleeve on the posterior end needle.
Length 25 mm
Diameter 22 G (0.72 mm)

Technical Specification :

Shelf life:
61 months.

Weight and Volume:
Estimated volume: 0.001 m3.
Estimated weight: 0.3 kg.



Laboratory Instrument India