



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Needle, vacuum tube, 22 G, box/100

**Product Code :**

LBNY-0001-21700015



**Description :**

Applicable for multi-sample draw  
Needle fits (screw in and out) needle holder  
Sterile, individually packed in a twist-apart plastic container  
Blood collection needle for use with vacuum tubes. For use in clinical phlebotomy.  
Material; Medical grade stainless steel  
Blood collection needle  
Pointed at both ends, with rubber sleeve on the posterior end needle.  
Length 25 mm  
Diameter 22 G (0.72 mm)

**Technical Specification :**

Shelf life:  
61 months.

Weight and Volume:  
Estimated volume: 0.001 m3.  
Estimated weight: 0.3 kg.



**Laboratory Instrument India**