



Email : sales@laboratoryinstrumentindia.com

Product Name :

Needle holder, vacuum tubes, box/100

Product Code :

LBNY-0001-21700016



Description :

Blood collection needle holder for use with vacuum tubes, and vacuum tube needles.
For use in clinical phlebotomy.
Material; clear polypropylene plastic
Fits all standard, 13mm, and 16mm diameter vacuum tubes, and all common needle units.
Disposable
Blood collection needle holder.

Technical Specification :

Shelf life:
36 months.

Weight and Volume:
Estimated volume: 0.0073 m3.
Estimated weight: 0.5 kg.



Laboratory Instrument India

